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

The image shows the 'Variable Editor for Liquid.DATA' dialog box. It has a 'Variables' list on the left with 'SNR' selected. A 'Repository' box on the right contains 'SNR'. Below these are 'Mandatory fields' including 'Identifier' (SNR), 'Name' (signal-noise ratio), 'Scope' (Peak), 'Type' (Real), 'Origin' (User Formula), 'Unit' (empty), and 'Formula' (HEIGHT/NOISE). There are also checkboxes for 'Display total', 'Mandatory', and 'Locked'. At the bottom are 'Optional fields' for 'Full Name', 'Full Unit', and 'Comment'. On the right side of the dialog are buttons for 'OK', 'Cancel', 'Help', 'New', 'Apply', 'Add', and 'Delete'.

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
Chromatogram properties...			
Main Information			
Identifier	Description	Value	Origin
NGROUP	#Groups	0	System : Number of groups in the chromatogram
NOISE	Noise [µV]	244	System : The noise is calculated between limits defined by "compute noise on" and "compute noise off" events
NOISE_ST	Noise SDEV	20.86	System : Estimation of channel noise standard deviation
NPEAKS	#Peaks	1	System : Number of integrated peaks
NPOINTS	#Pts	3901	System : Total number of acquired points of the chromatogram
OPERATOR	Operator Name	User1	System : Name of the user who performed the acquisition
PROCESS	Processing date	26/06/2011	System : Date of last processing

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