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EDINBURGH  
INSTRUMENTS

## The Supercontinuum Laser as a Flexible Source for Quasi-Steady State and Time Resolved Fluorescence Studies

Roger Fenske<sup>a,b</sup>, Dirk U. Näther<sup>a</sup>, Richard B. Dennisa, S. Desmond Smith<sup>a</sup>

<sup>a</sup>Edinburgh Instruments Ltd, 2 Bain Square, Kirkton Campus, Livingston, EH54 7DQ, UK;

<sup>b</sup>Heriot Watt University, Edinburgh, Scotland, EH14 4AS, UK

### ABSTRACT

Commercial Fluorescence Lifetime Spectrometers have long suffered from the lack of a simple, compact and relatively inexpensive broad spectral band light source that can be flexibly employed for both quasi-steady state and time resolved measurements (using Time Correlated Single Photon Counting [TCSPC]). This paper reports the integration of an optically pumped photonic crystal fibre, supercontinuum source<sup>1</sup> (Fianium model SC400PP) as a light source in Fluorescence Lifetime Spectrometers (Edinburgh Instruments FLS920 and Lifespec II), with single photon counting detectors (micro-channel plate photomultiplier and a near-infrared photomultiplier) covering the UV to NIR range. An innovative method of spectral selection of the supercontinuum source involving wedge interference filters is also discussed.

**Keywords:** Time Correlated Single Photon Counting, TCSPC, supercontinuum, fluorescence lifetimes

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TECHCOMP (CHINA) LTD.