Technical Note

**商品化 TCSPC 技术光谱仪的  
最新光源（摘要）**EDINBURGH  
INSTRUMENTS**New light sources for time correlated single photon counting in commercially available spectrometers**

Roger Fenske\*, Dirk U. Näther, Mark Goossens, S. Desmond Smith

Edinburgh Instruments, 2 Bain Square, Kirkton Campus, Livingston, UK, EH54 7DQ;

**ABSTRACT**

The method of Time Correlated Single Photon Counting requires high repetitive light sources ( $>100\text{kHz}$ ) with pulse widths of ideally less than approximately 20ps. While these light sources have been available for some time now in the form of Ti:Sapphire lasers, picosecond pulsed diode lasers ( $<90\text{ps}$ ) and light emitting diodes ( $<700\text{ps}$ ), they all have the drawback of either having no spectral tunability, or tunability over a very narrow spectral range (10nm-100nm).

While this is often sufficient for specific laboratory setups for measurements of fluorescence lifetimes, commercial Fluorescence Lifetime Spectrometers have suffered for a long time from the lack of the availability of simple, compact and relatively inexpensive broad spectral band light sources that can be employed for Time Correlated Single Photon Counting.

A new light source as an integral part of a commercial fluorescence lifetime spectrometer will be discussed that allows tunability over a wide spectral band of more than 500nm.

Keywords: Time Correlated Single Photon Counting, supercontinuum light source

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TECHCOMP (CHINA) LTD.